

CDK9 Rabbit mAb


ABZ2418

Product Overview


Name	CDK9 Rabbit mAb
Accession(Primary)	P50750
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:20-1:50; FC-1:50-1:100; IP-1:20

Target information (P50750)


Gene Name	CDK9
About Gene Symbol	



Western blot analysis of Cdk9 in K562, C6, 3T3, HeLa lysates using CDK9 antibody.



Immunocytochemistry analysis of Cdk9(green) in HeLa using Cdk9 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Cdk9 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.