

RASSF2 Rabbit mAb

ABZ2419


Product Overview

Name	RASSF2 Rabbit mAb
Accession(Primary)	P50749
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P50749)

Gene Name RASSF2

About Gene Symbol



Western blot analysis of RASSF2 in HeLa lysates using RASSF2 antibody.

