

VASP Rabbit mAb


ABZ2420

Product Overview

Name	VASP Rabbit mAb
Accession(Primary)	P50552
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 40 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:50-1:100

Target information (P50552)

Gene Name	VASP
About Gene Symbol	VASP (Vasodilator Stimulated Phosphoprotein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of VASP in HeLa lysates
using VASP antibody.