

Emerin Rabbit mAb


ABZ2421

Product Overview


Name	Emerin Rabbit mAb
Accession(Primary)	P50402
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 29 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; IF-1:20-1:50; FC-1:20-1:50

Target information (P50402)


Gene Name	EMD
About Gene Symbol	EMD (Emerin) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Emerin in 3T3, HeLa lysates using Emerin antibody.



Immunocytochemistry analysis of Emerin(green) in HeLa using Emerin antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Emerin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.