

MMP14 Rabbit mAb


ABZ2422

Product Overview


Name	MMP14 Rabbit mAb
Accession(Primary)	P50281
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 66 kDa; Observed MW: 66 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:100-1:200; IP-1:20

Target information (P50281)

Gene Name	MMP14
About Gene Symbol	MMP14 (Matrix Metalloproteinase 14) is a protease enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of MMP14 in CHO-K1 lysates using MMP14 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using MMP14 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.