

GNAQ Rabbit mAb


ABZ2423

Product Overview

Name	GNAQ Rabbit mAb
Accession(Primary)	P50148
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P50148)

Gene Name	GNAQ
About Gene Symbol	GNAQ (G Protein Subunit Alpha Q) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cel



Western blot analysis of GNAQ in rat Brain, C6, 3T3 lysates using GNAQ antibody.