

GNAQ Rabbit mAb


ABZ2424

Product Overview


Name	GNAQ Rabbit mAb
Accession(Primary)	P50148
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; FC-1:20-1:50; IP-1:20-1:50

Target information (P50148)

Gene Name	GNAQ
About Gene Symbol	GNAQ (G Protein Subunit Alpha Q) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cel



Western blot analysis of GNAQ in K562, rat Brain, C6, HeLa lysates using GNAQ antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using GNAQ antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.