

## CD39 Rabbit mAb

ABZ2425

### Product Overview

|                             |                                               |
|-----------------------------|-----------------------------------------------|
| <b>Name</b>                 | CD39 Rabbit mAb                               |
| <b>Accession(Primary)</b>   | P49961                                        |
| <b>Host Species</b>         | Rabbit                                        |
| <b>Clonality</b>            | Monoclonal                                    |
| <b>Isotype</b>              | IgG                                           |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                             |
| <b>Calculated MW</b>        | Calculated MW: 58 kDa; Observed MW: 34-70 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000                              |

### Target information (P49961)

**Gene Name** ENTPD1**About Gene Symbol**

