

Mre11 Rabbit mAb


ABZ2427

Product Overview


Name	Mre11 Rabbit mAb
Accession(Primary)	P49959
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 81 kDa; Observed MW: 81 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; FC-1:200-1:500

Target information (P49959)

Gene Name	MRE11
About Gene Symbol	



Western blot analysis of Mre11 in K562, rat Brain, C6, 3T3, Hela lysates using Mre11 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Mre11 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.