

Presenilin 2/AD5 Rabbit mAb


ABZ2433

Product Overview

Name	Presenilin 2/AD5 Rabbit mAb
Accession(Primary)	P49810
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:50-1:200; FC-1:20-1:50; IP-1:20

Target information (P49810)

Gene Name	PSEN2
About Gene Symbol	



Western blot analysis of Presenilin 2 in rat Brain, C6, 3T3, HeLa lysates using Presenilin 2 antibody.