

PLGF Rabbit mAb

ABZ2434


Product Overview

Name	PLGF Rabbit mAb
Accession(Primary)	P49763
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P49763)

Gene Name PGF

About Gene Symbol



Western blot analysis of PLGF in K562, C6, 3T3, HeLa lysates using PLGF antibody.

