

MCM2 Rabbit mAb


ABZ2437

Product Overview


Name	MCM2 Rabbit mAb
Accession(Primary)	P49736
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 125 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:200; ICC/IF-1:200-1:1000; FC-1:100-1:200

Target information (P49736)

Gene Name	MCM2
About Gene Symbol	



Western blot analysis of MCM2 in K562, C6, 3T3, HeLa lysates using MCM2 antibody.



Immunocytochemistry analysis of MCM2(green) in HeLa using MCM2 antibody, and DAPI(blue)