

MCM2 Rabbit mAb


ABZ2438

Product Overview


Name	MCM2 Rabbit mAb
Accession(Primary)	P49736
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 102 kDa; Observed MW: 125 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:20; ICC/IF-1:20-1:50; FC-1:20; IP-1:20

Target information (P49736)


Gene Name	MCM2
About Gene Symbol	



Western blot analysis of MCM2 in HeLa, A549, HL-60, U2OS lysates using MCM2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MCM2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of MCM2 (green)
in hela using MCM2 antibody, and DAPI(blue)