

HNF-4-alpha Rabbit mAb


ABZ2440

Product Overview

Name	HNF-4-alpha Rabbit mAb
Accession(Primary)	P49698
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:500-1:2000; FC-1:200-1:1000

Target information (P49698)

Gene Name	Hnf4A
About Gene Symbol	



Western blot analysis of HNF4alpha in Raw264.7
lysates using HNF 4 alpha antibody.