

HNF-4-alpha Rabbit mAb


ABZ2441

Product Overview


Name	HNF-4-alpha Rabbit mAb
Accession(Primary)	P49698
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20; ChIP-1:20; IHC-1:50-1:100

Target information (P49698)

Gene Name	Hnf4A
About Gene Symbol	



Western blot analysis of HNF4alpha in mouse liver lysates using HNF 4 alpha antibody.



Western blot analysis of HNF4alpha in K562, C6, HeLa lysates using HNF4alpha antibody.