

ALDH7A1 Rabbit mAb


ABZ2443

Product Overview


Name	ALDH7A1 Rabbit mAb
Accession(Primary)	P49419
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:50-1:100

Target information (P49419)


Gene Name	ALDH7A1
About Gene Symbol	



Western blot analysis of ALDH7A1 in mouse liver lysates using ALDH7A1 antibody.



Western blot analysis of ALDH7A1 in HeLa lysates using ALDH7A1 antibody.



Immunocytochemistry analysis of
ALDH7A1(green) in HeLa using ALDH7A1
antibody, and DAPI(blue)