

Fatty Acid Synthase Rabbit mAb


ABZ2446

Product Overview

Name	Fatty Acid Synthase Rabbit mAb
Accession(Primary)	P49327
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 273 kDa; Observed MW: 273 kDa
Recommended Dilution	WB-1:1000-1:10000; IHC-1:100-1:10000; ICC/IF-1:200-1:500; FC-1:20-1:50; IP-1:20-1:50


Target information (P49327)

Gene Name	FASN
About Gene Symbol	FASN (Fatty Acid Synthase) is a metabolic dehydrogenase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of Fatty Acid Synthase in 293T, A549 lysates using Fatty Acid Synthase antibody.



Western blot analysis of Fatty Acid Synthase in HeLa lysates using Fatty Acid Synthase antibody.