

LMAN1 Rabbit mAb


ABZ2447

Product Overview


Name	LMAN1 Rabbit mAb
Accession(Primary)	P49257
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100; FC-1:20

Target information (P49257)

Gene Name	LMAN1
About Gene Symbol	LMAN1 (Lectin, Mannose Binding 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of LMAN1 in HeLa, A549, HL-60, U251, U87-MG lysates using LMAN1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using LMAN1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.