

MK2 Kinase 2 Rabbit mAb


ABZ2448

Product Overview


Name	MK2 Kinase 2 Rabbit mAb
Accession(Primary)	P49137
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P49137)

Gene Name	MAPKAPK2
About Gene Symbol	MAPKAPK2 (MAPK Activated Protein Kinase 2) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of MK2 in mouse heart, mouse skeletal muscle, rat skeletal muscle lysates using MAPKAP Kinase 2 antibody.



Western blot analysis of MK2 in HL-60, U2OS lysates using MK2 antibody.