

Paxillin Rabbit mAb


ABZ2450

Product Overview


Name	Paxillin Rabbit mAb
Accession(Primary)	P49023
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 65 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:50-1:100; FC-1:20; IP-1:20

Target information (P49023)


Gene Name	PXN
About Gene Symbol	PXN (Paxillin) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of Paxillin in mouse lung lysates using Paxillin antibody.



Western blot analysis of Paxillin in C6 lysates using Paxillin antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Paxillin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Paxillin (green) in HeLa using Paxillin antibody, and DAPI (blue).