


PI3 Kinase p110 gamma Rabbit mAb

ABZ2451

Product Overview

Name	PI3 Kinase p110 gamma Rabbit mAb
Accession(Primary)	P48736
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 126 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P48736)

Gene Name PIK3CG**About Gene Symbol**

Western blot analysis of PI3 Kinase p110 gamma in K562 lysates using PI3 Kinase p110 gamma antibody.

