

SOX9 Rabbit mAb


ABZ2458

Product Overview


Name	SOX9 Rabbit mAb
Accession(Primary)	P48436
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:2000; IHC:1:100-1:200; ICC/IF-1:20-1:50; FC-1:20; IP-1:20

Target information (P48436)

Gene Name	SOX9
About Gene Symbol	SOX9 (SRY-Box Transcription Factor 9) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of SOX9 in 3T3, HeLa lysates using SOX9 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using SOX9 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.