

eIF1A Rabbit mAb


ABZ2461

Product Overview


Name	eIF1A Rabbit mAb
Accession(Primary)	P47813
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100; IP-1:10-1:100

Target information (P47813)


Gene Name	EIF1AX
About Gene Symbol	




Immunocytochemistry analysis of eIF1A (green) in 293 using eIF1A antibody, and DAPI (blue).



Western blot analysis of eIF1A in mouse brain, rat brain lysates using eIF1A antibody.



Western blot analysis of eIF1A in LnCap, K562, U251MG, HepG2, MCF-7, HeLa lysates using eIF1A antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using eIF1A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.