

RAP1GAP Rabbit mAb


ABZ2463

Product Overview


Name	RAP1GAP Rabbit mAb
Accession(Primary)	P47736
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:20-1:100

Target information (P47736)


Gene Name	RAP1GAP
About Gene Symbol	




Western blot analysis of RAP1GAP in mouse brain lysates using RAP1GAP antibody.



Western blot analysis of RAP1GAP in rat Brain, HeLa lysates using RAP1GAP antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using RAP1GAP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of RAP1GAP(green) in Hela using RAP1GAP antibody, and DAPI(blue).