

RAP1GAP Rabbit mAb


ABZ2464

Product Overview


Name	RAP1GAP Rabbit mAb
Accession(Primary)	P47736
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:20-1:50; FC-1:20; IP-1:20

Target information (P47736)


Gene Name	RAP1GAP
About Gene Symbol	



Western blot analysis of RAP1GAP in A549 lysates using RAP1GAP antibody.



Western blot analysis of RAP1GAP in C6, 3T3 lysates using RAP1GAP antibody.



Immunocytochemistry analysis of
RAP1GAP(green) in HeLa using RAP1GAP
antibody, and DAPI(blue)