

YAP1 Rabbit mAb


ABZ2466

Product Overview


Name	YAP1 Rabbit mAb
Accession(Primary)	P46937
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 54 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:100-1:500; FC-1:20; IP-1:20

Target information (P46937)

Gene Name	YAP1
About Gene Symbol	YAP1 (Yes1 Associated Transcriptional Regulator) is a transcriptional coactivator involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of YAP1 in HeLa lysates using YAP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using YAP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.