


NEDD4 Rabbit mAb

ABZ2467


Product Overview

Name	NEDD4 Rabbit mAb
Accession(Primary)	P46934
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 149 kDa; Observed MW: 115 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:20-1:50; IP-1:20-1:50

Target information (P46934)

Gene Name NEDD4**About Gene Symbol**

Western blot analysis of NEDD4 in K562, HCT116, 3T3, C6 lysates using NEDD4 antibody.



Western blot analysis of NEDD4 in HeLa, A549, HL-60, U2OS, U87-MG lysates using NEDD4 antibody.