


BRCC36 Rabbit mAb

ABZ2468

Product Overview

Name	BRCC36 Rabbit mAb
Accession(Primary)	P46736
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:500

Target information (P46736)

Gene Name BRCC3**About Gene Symbol**

Western blot analysis of BRCC36 in rat Brain, C6, 3T3, HeLa lysates using BRCC36 antibody.

