

MEK3/MEK6 Rabbit mAb


ABZ2470

Product Overview


Name	MEK3/MEK6 Rabbit mAb
Accession(Primary)	P46734
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 39,37 kDa; Observed MW: 39,37 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; IF-1:200-1:1000; FC-1:100-1:200; IP-1:20-1:100

Target information (P46734)


Gene Name	MAP2K3/MAP2K6
About Gene Symbol	MAP2K3 (Mitogen-Activated Protein Kinase Kinase 3) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of MEK3/MEK6 (green) in HCT116 using MEK3/MEK6 antibody, and DAPI (blue).



Western blot analysis of MEK3/MEK6 in 3T3 lysates using MEK3/MEK6 antibody. Observed band size: 39, 37 kDa.



Western blot analysis of MEK3/MEK6 in HeLa, A549, HL-60, U2OS lysates using MEK3/MEK6 antibody