

RYP Rabbit mAb


ABZ2471

Product Overview


Name	RYP Rabbit mAb
Accession(Primary)	P34925
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 68 kDa; Observed MW: 72 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:20-1:100; IP-1:10-1:100

Target information (P34925)


Gene Name	RYP
About Gene Symbol	



Western blot analysis of RYK in Hela, A549, HL-60, U2OS, U87-MG lysates using RYK antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using RYK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of RYK (green) in
hela using RYK antibody, and DAPI(blue)