

p27 Kip 1 Rabbit mAb


ABZ2473

Product Overview


Name	p27 Kip 1 Rabbit mAb
Accession(Primary)	P46414
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 22 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:1000-1:5000; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (P46414)


Gene Name	Cdkn1B
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse lung using p27 Kip 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of p27 KIP 1 in 3T3 lysates using p27 KIP 1 antibody.



Western blot analysis of p27 KIP 1 in HeLa, A549, HL-60 lysates using p27 KIP 1 antibody