

## BAT3 Rabbit mAb


ABZ2474

### Product Overview


|                             |                                                                         |
|-----------------------------|-------------------------------------------------------------------------|
| <b>Name</b>                 | BAT3 Rabbit mAb                                                         |
| <b>Accession(Primary)</b>   | P46379                                                                  |
| <b>Host Species</b>         | Rabbit                                                                  |
| <b>Clonality</b>            | Monoclonal                                                              |
| <b>Isotype</b>              | IgG                                                                     |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                       |
| <b>Calculated MW</b>        | Calculated MW: 119 kDa; Observed MW: 150 kDa                            |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:50-1:100;<br>ICC/IF-1:200-1:500; FC-1:20-1:50 |

### Target information (P46379)

|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | BAG6                                                                                                                                                                                                                                                            |
| <b>About Gene Symbol</b> | BAG6 (BAG Cochaperone 6) is a molecular chaperone involved in apoptosis. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of BAT3/BAG6 in K562, rat Brain, C6, 3T3 lysates using BAT3 antibody.



Immunocytochemistry analysis of BAT3/BAG6 (green) in MCF-7 using BAT3/BAG6 antibody, and DAPI (blue)