

CrkL Rabbit mAb


ABZ2475

Product Overview



Name	CrkL Rabbit mAb
Accession(Primary)	P46109
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 34 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100

Target information (P46109)

Gene Name	CRKL
About Gene Symbol	CRKL (CRK Like Proto-Oncogene, Adaptor Protein) is a signaling adaptor protein involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of CrkL (green) in K562 using CrkL antibody, and DAPI (blue).



Western blot analysis of CrkL in Jurkat lysates using CrkL antibody.