

RanGAP1 Rabbit mAb


ABZ2476

Product Overview


Name	RanGAP1 Rabbit mAb
Accession(Primary)	P46060
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 64 kDa; Observed MW: 68,82 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:20

Target information (P46060)


Gene Name	RANGAP1
About Gene Symbol	




Western blot analysis of RanGAP1 in 3T3 lysates using RanGAP1 antibody.



Western blot analysis of RanGAP1 in HeLa, A549, HL-60, U2OS, C6 lysates using RanGAP1 antibody.



Immunocytochemistry analysis of RanGAP1 (green) in HeLa using RanGAP1 antibody, and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using RanGAP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.