

JNK2 Rabbit mAb


ABZ2477

Product Overview


Name	JNK2 Rabbit mAb
Accession(Primary)	P45984
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:2000; IHC:1:20-1:50; FC-1:20; IP-1:20

Target information (P45984)

Gene Name	MAPK9
About Gene Symbol	



Western blot analysis of JNK2 in C6, 3T3 lysates using JNK2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using JNK2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.