

## JNK1+JNK2+JNK3 Rabbit mAb


ABZ2479

### Product Overview


|                             |                                                               |
|-----------------------------|---------------------------------------------------------------|
| <b>Name</b>                 | JNK1+JNK2+JNK3 Rabbit mAb                                     |
| <b>Accession(Primary)</b>   | P45983                                                        |
| <b>Host Species</b>         | Rabbit                                                        |
| <b>Clonality</b>            | Monoclonal                                                    |
| <b>Isotype</b>              | IgG                                                           |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                             |
| <b>Calculated MW</b>        | Calculated MW: 48 kDa; Observed MW: 46,54 kDa                 |
| <b>Recommended Dilution</b> | WB-1:20000-1:100000; IF: 1:50; FC-1:50-1:100;<br>IP-1:20-1:50 |

### Target information (P45983)

|                          |                                                                                                                                                                                                                                                                                                         |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | MAPK8                                                                                                                                                                                                                                                                                                   |
| <b>About Gene Symbol</b> | MAPK8 (Mitogen-Activated Protein Kinase 8) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of JNK1/2/3 in mouse brain lysates using JNK antibody.



Western blot analysis of JNK1/2/3 in Jurkat, rat Brain, C6, CHO-K1, Hela lysates using JNK1/2/3 antibody.