

JNK1+JNK2+JNK3 Rabbit mAb


ABZ2480

Product Overview


Name	JNK1+JNK2+JNK3 Rabbit mAb
Accession(Primary)	P45983
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 46,54 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100-1:200; FC-1:100-1:200; IP-1; 20-1:50

Target information (P45983)

Gene Name	MAPK8
About Gene Symbol	MAPK8 (Mitogen-Activated Protein Kinase 8) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of JNK1/2/3 in C6, 3T3 lysates using JNK antibody.



Immunocytochemistry analysis of JNK1/2/3 (green) in HeLa using JNK1/2/3 antibody, and DAPI (blue)