


ABCG1 Rabbit mAb

ABZ2482

Product Overview

Name	ABCG1 Rabbit mAb
Accession(Primary)	P45844
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 76 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100

Target information (P45844)

Gene Name ABCG1**About Gene Symbol**

Western blot analysis of ABCG1 in HeLa lysates using ABCG1 antibody.

