

Nicotinic Acetylcholine Receptor alpha 4 Rabbit mAb


ABZ2483

Product Overview


Name	Nicotinic Acetylcholine Receptor alpha 4 Rabbit mAb
Accession(Primary)	P43681
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (P43681)

Gene Name	CHRNA4
About Gene Symbol	



Western blot analysis of Nicotinic Acetylcholine Receptor alpha 4/CHRNA4 in 3T3 lysates using Nicotinic Acetylcholine Receptor alpha 4 antibody.



Western blot analysis of Nicotinic Acetylcholine Receptor alpha 4/CHRNA4 in HeLa, CHO-K1, C6, Ramos lysates using Nicotinic Acetylcholine Receptor alpha 4/CHRNA4 antibody.