

ZAP70 Rabbit mAb

ABZ2484


Product Overview

Name	ZAP70 Rabbit mAb
Accession(Primary)	P43403
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:100; FC-1:20-1:50

Target information (P43403)

Gene Name ZAP70

About Gene Symbol



Western blot analysis of ZAP70 in Jurkat lysates using ZAP70 antibody.

