


GPD2 Rabbit mAb

ABZ2485


Product Overview

Name	GPD2 Rabbit mAb
Accession(Primary)	P43304
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 81 kDa; Observed MW: 81 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P43304)

Gene Name GPD2**About Gene Symbol**

Western blot analysis of GPD2 in rat testis, rat skeletal muscle lysates using GPD2 antibody.



Western blot analysis of GPD2 in K562, HeLa
lysates using GPD2 antibody.