


MSH2 Rabbit mAb

ABZ2486

Product Overview

Name	MSH2 Rabbit mAb
Accession(Primary)	P43246
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 105 kDa; Observed MW: 105 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P43246)

Gene Name MSH2**About Gene Symbol**

Western blot analysis of MSH2 in rat Brain, CHO-K1, HeLa lysates using MSH2 antibody.

