

LIS1 Rabbit mAb


ABZ2489

Product Overview

Name	LIS1 Rabbit mAb
Accession(Primary)	P43034
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P43034)

Gene Name	PAFAH1B1
About Gene Symbol	PAFAH1B1 (Platelet Activating Factor Acetylhydrolase 1b Regulatory Subunit 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to el



Western blot analysis of LIS1 in C6, CHO-K1
lysates using LIS1 antibody.