

p18 INK4c/CDKN2C Rabbit mAb


ABZ2490

Product Overview


Name	p18 INK4c/CDKN2C Rabbit mAb
Accession(Primary)	P42773
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 18 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (P42773)


Gene Name	CDKN2C
About Gene Symbol	CDKN2C (Cyclin Dependent Kinase Inhibitor 2C) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of p18 INK4c/CDKN2C in 3T3 lysates using CDKN2C antibody.



Western blot analysis of p18 INK4c/CDKN2C in C6 lysates using p18 INK4c/CDKN2C antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using p18 INK4c antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.