

CDKN2A/p16INK4a Rabbit mAb


ABZ2491

Product Overview


Name	CDKN2A/p16INK4a Rabbit mAb
Accession(Primary)	P42771
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:20-1:100; FC-1:20

Target information (P42771)


Gene Name	CDKN2A
About Gene Symbol	CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) is a tumor suppressor protein involved in cell cycle regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CDKN2A/p16INK4a in 293 lysates using CDKN2A antibody.



Western blot analysis of CDKN2A/p16INK4a in 3T3, HeLa lysates using CDKN2A/p16INK4a antibody.



Immunocytochemistry analysis of
CDKN2A/p16INK4a(green) in HeLa using
CDKN2A/p16INK4a antibody, and DAPI(blue)