

Tec Rabbit mAb


ABZ2494

Product Overview


Name	Tec Rabbit mAb
Accession(Primary)	P42680
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 74 kDa; Observed MW: 58-75 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:500; FC-1:20-1:100

Target information (P42680)

Gene Name	TEC
About Gene Symbol	



Western blot analysis of Tec in K562 lysates using Tec antibody.



Immunocytochemistry analysis of Tec(green) in Hela using Tec antibody, and DAPI(blue)