

Caspase 3 Rabbit mAb


ABZ2497

Product Overview


Name	Caspase 3 Rabbit mAb
Accession(Primary)	P42574
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 32 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; IP-1:50-1:100

Target information (P42574)

Gene Name	CASP3
About Gene Symbol	CASP3 (Caspase 3) is a cysteine protease (caspase) involved in apoptosis. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Caspase3 in K562, rat Brain, C6, 3T3, HeLa lysates using Caspase 3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Caspase3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.