

PI3 Kinase p110 alpha Rabbit mAb


ABZ2498

Product Overview

Name	PI3 Kinase p110 alpha Rabbit mAb
Accession(Primary)	P42336
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 124 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:2000; IF: 1:50; FC-1:20-1:50; IP-1:20-1:50


Target information (P42336)

Gene Name	PIK3CA
About Gene Symbol	



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of PI3-Kinase catalytic subunit alpha in Jurkat, MCF-7, Raw264.7 lysates using PI3-Kinase catalytic subunit alpha antibody.



Western blot analysis of PI3 Kinase p110 alpha in Jurkat lysates using PI3 Kinase p110 alpha antibody.