

Glutamate receptor 2 Rabbit mAb


ABZ2499

Product Overview


Name	Glutamate receptor 2 Rabbit mAb
Accession(Primary)	P42262
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 99 kDa; Observed MW: 99 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:50-1:100; IHC-1:50-1:100

Target information (P42262)

Gene Name	GRIA2
About Gene Symbol	



Western blot analysis of Ionotropic Glutamate receptor 2 in C6 lysates using Ionotropic Glutamate Receptor 2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Ionotropic Glutamate receptor 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.