

STAT5 Rabbit mAb


ABZ2500

Product Overview


Name	STAT5 Rabbit mAb
Accession(Primary)	P42230
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 91 kDa; Observed MW: 91 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:50-1:100

Target information (P42230)

Gene Name	Stat5A
About Gene Symbol	



Western blot analysis of STAT5a/b in K562, C6, 3T3, HeLa lysates using STAT5 antibody.



Immunocytochemistry analysis of STAT5A/B(green) in HeLa using STAT5A/B antibody, and DAPI(blue)